

# microFET 6™

## microFET 6™ Range of Motion Tester

- Can be used as either a Dual or Single inclinometer
- Ergonomically designed to be lightweight and easy to use
- Designed for optimal instrument placement and spinal clearance
- Easy to read LCD screen
- Precision casing for durability and extended instrument life
- High reliability, durability and low incidence of repairs – for low total cost
- For use with AMA or DRE testing protocols
- For use as standalone or with our clinical software



EVALUATION TOOLS TO MEASURE, OBJECTIFY AND DOCUMENT HUMAN PERFORMANCE AND CARE

## microFET FAMILY



Manual Muscle Tester



MT/ ROM



Pinch & Grip



Intelligent ROM

# microFET 6™



## microFET 6™ Specifications

- Size: Main unit: 2.5" x 4" Secondary unit: 2.5" x 2"
- Weight: Main unit: .12 oz secondary unit: 4 oz.
- 1° increments from +/- 0° to +/- 180°
- Controls: Reset & Enter
- Accurate to within 1°
- 1 year warranty (Extended Warranties Available)
- Self-activating "sleep" mode to extend battery life

## microFET 6 Product Includes:

- microFET6 - Primary and Secondary Devices
- Connecting Cord
- User Manual
- Protocol Manual
- Carrying Case
- ROM Wallchart
- 3 test recording pads for manual testing
- Calibration Certificate